

YOUR CONTACT
HERE

YOUR LOGO HERE



- 3 PLY EAR-LOOP
- WITHOUT GLASS FIBERS
- HYPOLLERGENIC
- LOW RESISTANCE TO BREATHING
- NOSE BAR ADAPTABLE
- HIGH FILTRATION CAPACITY
- PERFECT FITTING

50 pcs/box 40 boxes/case



Disposable Protective Mask
LISTING # 10063057
PRODUCT CODE: QKR

GREAT QUALITY MASK
YOUR BRAND
YOUR LOGO



U.S. BASED
CUSTOMER SERVICE



US BASED
GRAPHIC DESIGN



TURN AROUND YOUR DESIGN
IN NO TIME

ADULT MASK TESTING REPORT

№ 轻委2020-06-0730

LC2006094-2

广州质量监督检测研究院
Guangzhou Quality Supervision and Testing Institute (GQT)

Page 2 of 3

检验检测报告 Test Report

第 2 页 共 3 页

序号 No.	检验项目 Test Items	单位 Unit	标准要求 Reference Requirement	检验结果 Test Results	单项评价 Assess
			---	---	
1	通气阻力 Ventilation resistance	Pa	儿童/老人: ≤50 其他: ≤70 Children / The olders: ≤50 Others: ≤70	1#:18.6 2#:17.8 3#:17.0	合格 Pass
2	颗粒物过滤效率 Particle filtration efficiency	%	A级: ≥50 Level A: ≥50	1#:96.1 2#:97.4 3#:96.7	合格 Pass
3	细菌过滤效率 (BFE) Bacterial filtration efficiency (BFE)	%	口罩的细菌过滤效率应不小于95 The filtration efficiency of the mask should be no less than 95	1#: 99.7 2#: 99.5 3#: 99.9	合格 Pass

KID'S MASK TESTING REPORT



Study Number 1323308-S01
Bacterial Filtration Efficiency (BFE)
and Differential Pressure (Delta P) Final Report

Results:

Test Article Number	Percent BFE (%)
1	>99.9
2	>99.9
3	>99.9
4	99.7
5	>99.9

Test Article Number	Delta P (mm H ₂ O/cm ³)	Delta P (Pa/cm ³)
1	4.1	40.6
2	4.1	40.1
3	4.2	41.5
4	4.4	43.0
5	4.4	43.1

The filtration efficiency percentages were calculated using the following equation:

$$\% BFE = \frac{C - T}{C} \times 100$$

C = Positive control average

T = Plate count total recovered downstream of the test article

Note: The plate count total is available upon request

Establishment Registration & Device Listing

[FDA Home](#) |
 [Medical Devices](#) |
 [Databases](#)

[New Search](#)

[Back To Search Results](#)

Proprietary Name:	Fushentang; Kombang
Classification Name:	FACE MASK (EXCEPT N95 RESPIRATOR) FOR GENERAL PUBLIC/HEALTHCARE PERSONNEL PER IIE GUIDANCE
Product Code:	QKR
Device Class:	Not Classified
Registered Establishment Name:	GUANGZHOU KANGBANG MEDICAL EQUIPMENT CO., LTD.
Registered Establishment Number:	3016640643
Owner/Operator:	Guangzhou Kangbang Medical Equipment Co., Ltd.
Owner/Operator Number:	10063057
Establishment Operations:	Manufacturer